Titanized Mesh Implants

TiLOOP® PRO A

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
TiLOOP® PRO A

Description

Non-absorbable prosthetic mesh of monofilament polypropylene thread with surrounding titanized coating. TiLOOP® mesh implants are available in various basis weights with defined and customized shapes and dimensions.

TiLOOP® mesh implants were developed for the treatment of urology, gynecology and surgical patients with various conditions that have to date been treated with surgery involving mesh implants.

TiLOOP® mesh implants are used in particular for supporting and strengthening connective tissue structures and ligaments.

These Instructions for Use are also valid for other TiLOOP® PRO A product variants.

Indication

The TiLOOP® PRO pelvic floor reconstruction systems are available as a titanized mesh implants for strengthening and stabilising the fascia structure of the pelvic floor during surgical procedures. The various cuts cover the surgical procedure for rectocele, enterocele, and cystocele as well as the treatment of vaginal vault prolapse.

Mode of action

TiLOOP® mesh implants are used to strengthen tissue or bridge connective tissue structures and ligaments. The deposition of a thin fibrous tissue layer, whose cells grow through the pores of the mesh fabric, results in the complete enclosure of the mesh material by the body's own tissue.

Experimental investigations have shown that the titanium mesh material of TiLOOP® mesh implants remains soft and flexible.

The titanium coating of the synthetic material will not affect radiological diagnostic techniques (X-ray, CT, MRI and ultrasonic examinations).

Warnings

- ⚠ The physician should be experienced with the surgical procedures and techniques.
- ⚠ Vaginal or urinary infections should be treated before implantation.
- ⚠ Contaminated or infected wounds should be treated according to good surgical practice.
- ⚠ Risks and benefits of an implantation with known systemic steroid therapy must be carefully considered.
- ⚠ Risks and benefits of an implantation with known blood clotting dysfunction must be carefully considered.
- ⚠ Take care to avoid perforation of the urethra, bladder, bowel or vessel or damage of nerves during use of implantation instruments.
- ⚠ In case of injury immediate surgical intervention is necessary.
- ⚠ To exclude bladder perforation after mesh implantation, cystoscopy should be performed.
- ⚠ Patients have to be observed with respect to signs of secondary hemorrhage.
- ⚠ The patients will be examined for ureteral lesion/stenosis after implantation using suitable methods (for example ultrasound).
- ⚠ TiLOOP® mesh consists of polypropylene mesh that should not be damaged. When using instruments pay attention not to damage the mesh.
- ⚠ Do not use the product if the packaging is open or damaged or should you have doubts about the sterility or if the use by date has expired.
- ⚠ Each product is packed separately and is supplied EO-sterilised and pyrogen-free. It is only to be used once. Do not reuse, reprocess or resterilise. A reuse, reprocessing or resterilisation of disposable products can result in a reduced performance or a loss of their functionality. A reuse of disposable products can result in exposure to pathogens such as viruses, bacteria, fungi or prions.

Recommended implantation method using the TiLOOP® Total Application Set

In surgery, the patient will be prepared under strictly sterile conditions in accordance with clinical standards for the surgery.

Step 01: patient preparation

- The repair of the prolapse using the TiLOOP® pelvic floor reconstruction mesh will be carried out under regional anaesthesia or general anaesthesia.
- In surgery, the patient will be placed in the lithotomy position for surgery; the legs will be bent at the hips and placed on supports, and the buttocks will be moved to the edge of the surgical table or positioned slightly past the edge.
- Padding should be placed under the renal pelvis area so as to prevent kinking of the ureter or the ureter from being positioned in the operative region. Alternatively or additionally, it is recommended to use a temporary ureteral stent in step 3.
- The bladder should be emptied prior to the intervention. The use of a bladder catheter is recommended. The physician may decide whether a suprapubic or transurethral catheter should be used.

Step 02: preparation of the TiLOOP® mesh

- Ensure that the sterile barrier is intact before removing the TiLOOP® mesh from the packaging.
- Then open the packaging of the TiLOOP® mesh in a sterile environment.
- Verify that the contents of the packaging are complete and undamaged. Use another package if components are missing or show damage.

Step 03: dissection of the vagina, endopelvic fascia, and ischial spine

- Perform an anterior medial colpotomy until about 3 cm away from the external urinary meatus.
- The bladder will be separated by dissection; at this point it is important to ensure that blood flow to the vagina remains uninterrupted by performing full-thickness dissection.
- The dissection will be adequate when the exposure has been made laterally up to the vaginal sulcus.
- Then identify and expose the ischial spine.
- Repeat this process on the opposite side under careful interruption of blood flow.
- Starting from the ischial spine, dissect the sacrospinous ligament in the medial and dorsal direction.
- Repeat this process on the opposite side under careful interruption of blood flow.
- Alternatively or additionally to the padded position of the patient, it is recommended to use a temporary ureteral stent, which can be removed after fixing the mesh arms dorsally.

Step 04: preparation for instrument passage

- Mark the sites for the skin incisions for the passage of the TiLOOP® Instruments PLUS.
- The first skin incision should be made at the genitofemoral fold on the lateral margin of the descending pubic ramus at the level of the clitoris underneath the attachment of the adductor longus muscle.
- The second incision will be made 1-2 cm laterally and 2-3 cm caudally from the first.
- Carry out this process on the opposite side as well.
- A third incision will be made 3 cm lateral and dorsal to the anus.
- Carry out this process on the opposite side as well.

Step 05: inserting the instruments and TiLOOP® mesh arms

- Introduce the TiLOOP® Total Catheter over the TiLOOP® Instruments PLUS and insert it through the anterior incision. It will closely circumscribe the ascending pubic ramus. Hold the tip of the instrument in your fingers so as to prevent injuring the bladder.
- The TiLOOP® Total Catheter will remain in the body temporarily after the TiLOOP® Instruments PLUS has been removed.
- Insert the TiLOOP® Total Snare through the TiLOOP® Total Catheter using the free loop until the free loop exits the TiLOOP® Total Catheter distal to the vaginal incision. The TiLOOP® Total Snare will be left as is until all fixation points are demonstrated and prepared.
- Repeat this process on the opposite side.
- Introduce another TiLOOP® Total Catheter over the TiLOOP® Instruments PLUS; inserted it into the second incision and advance it past the obturator foramen into the dorsocaudal corner.
- Continue to advance the instrument at the pubic bone along the ischial spine behind the obturator internus muscle until the tip of the needle emerges laterocranial to the ischial spine. Insert another TiLOOP® Total Snare through the TiLOOP® Total Catheter using the free loop until the free loop exits the TiLOOP® Total Catheter distal to the vaginal incision. The TiLOOP® Total Snare will also be left as is.
- Repeat this process on the opposite side.
- Introduce another TiLOOP® Total Catheter over the TiLOOP® Instruments PLUS and insert it 3 cm lateral and dorsal to the anus. Then guide it through the ischioanal fossa, and pass it through the sacrospinous ligament dorsally at least 2.5-3 cm medial to the spine. Insert another TiLOOP® Total Snare through the TiLOOP® Total Catheter using the free loop until the free loop exits the TiLOOP® Total Catheter distal to the vaginal incision. The TiLOOP® Total Snare will also be left as is.
- Repeat this process on the opposite side.
- The TiLOOP® PRO A mesh can be grasped at the free loop by turning over a mesh arm. Next extract the mesh arm held by the TiLOOP® Total Snare from the body through the TiLOOP® Total Catheter using the TiLOOP® Total Snare.
- Repeat the steps described above for each TiLOOP® Total Catheter and each TiLOOP® Total Snare.
- Remove the TiLOOP® Total Catheters from the patient's tissue once the mesh is in the correct location and has been adjusted.

Step 06: securing the TiLOOP® mesh

- Ensure that the mesh is spread out flat and is without tension.
- The mesh will be fixed in this position with an absorbable suture at the vesical-urethral junction.
- In cases where the uterus is intact, the craniodorsal section of the mesh will be fixed to the cervix. After prior hysterectomy or simultaneous hysterectomy performed with this surgery, the mesh will be attached to both sides of the vaginal apex.

Step 07: suture

- Last, the vaginal wound will be fully closed and sutured using absorbable sutures.
- Tamponade the vagina after the surgery.

Recommended implantation method without using the TiLOOP® Total Application Set

In surgery, the patient will be prepared under strictly sterile conditions in accordance with clinical standards for the surgery.

Step 01: patient preparation

- The repair of the prolapse using the TiLOOP® pelvic floor reconstruction mesh will be carried out under regional anaesthesia or general anaesthesia.
- In surgery, the patient will be placed in the lithotomy position for surgery; the legs will be bent at the hips and placed on supports, and the buttocks will be moved to the edge of the surgical table or positioned slightly past the edge.
- Padding should be placed under the renal pelvis area so as to prevent kinking of the ureter or the ureter from being positioned in the operative region. Alternatively or additionally, it is recommended to use a temporary ureteral stent in step 3.
- The bladder should be emptied prior to the intervention. The use of a bladder catheter is recommended. The physician may decide whether a suprapubic or transurethral catheter should be used.

Step 02: preparation of the TiLOOP® mesh

- Ensure that the sterile barrier is intact before removing the TiLOOP® mesh from the packaging.
- Then open the packaging of the TiLOOP® mesh in a sterile environment.
- Verify that the contents of the packaging are complete and undamaged. Use another package if components are missing or show damage.

Step 03: dissection of the vagina, endopelvic fascia, and ischial spine

- Perform an anterior medial colpotomy until about 3 cm away from the external urinary meatus.
- The bladder will be separated by dissection; at this point it is important to ensure that blood flow to the vagina remains uninterrupted by performing full-thickness dissection.
- The dissection will be adequate when the exposure has been made laterally up to the vaginal sulcus.
- Then identify and expose the ischial spine.
- Repeat this process on the opposite side under careful interruption of blood flow.
- Starting from the ischial spine, dissect the sacrospinous ligament in the medial and dorsal direction.
- Repeat this process on the opposite side under careful interruption of blood flow.
- Alternatively or additionally to the padded position of the patient, it is recommended to use a temporary ureteral stent, which can be removed after fixing the mesh arms dorsally.

Step 04: preparation for instrument passage

- Mark the sites for the skin incision for the passage of the instruments (needle applicators, tunneler).
- The first skin incision should be made at the genitofemoral fold on the lateral margin of the descending pubic ramus at the level of the clitoris underneath the attachment of the adductor longus muscle.
- The second incision will be made 1-2 cm laterally and 2-3 cm caudally from the first.
- Carry out this process on the opposite side as well.
- A third incision will be made 3 cm lateral and dorsal to the anus.
- Carry out this process on the opposite side as well.

Step 05: inserting the instruments and TiLOOP® mesh arms

- Insert the helical tunneler through the anterior incision. It will closely circumscribe the ascending pubic ramus. Hold the tip of the instrument in your fingers so as to prevent injuring the bladder.
- Hook the anterior arm of the TiLOOP® PRO A mesh into the eye of the helical tunneler (the mechanism for attaching the mesh may vary depending on the instrument used), and carefully extract the mesh arm.
- Repeat this process on the opposite side.
- Insert the curved tunneler through the second incision and advance it past the obturator foramen in the dorsocaudal corner.
- Continue to advance the instrument at the pubic bone along the ischial spine behind the obturator internus muscle until the tip of the needle emerges laterocranial to the ischial spine.
- Now hook the middle arm of the TiLOOP® PRO A mesh into the eye of the tunneler (the mechanism for attaching the mesh may vary depending on the instrument used), and carefully extract the mesh arm.
- Repeat this process on the opposite side.
- Insert the curved tunneler 3 cm lateral and dorsal to the anus. From there, guide it through the ischioanal fossa, and pass it through the sacrospinous ligament dorsally, 2.5-3 cm medial to the spine.
- Hook the dorsal arm of the TiLOOP® PRO A mesh into the eye of the tunneler (the mechanism for attaching the mesh may vary depending on the instrument used), and carefully extract the mesh arm.
- Repeat this process on the opposite side.
- It may be indicated to use non-absorbable sutures to guide the mesh arms instead of hooking into the mesh arms directly and guiding them using the hooked instrument.

Step 06: securing the TiLOOP® mesh

- Ensure that the mesh is spread out flat and is without tension.
- The mesh will be fixed in this position with an absorbable suture at the vesical-urethral junction.
- In cases where the uterus is intact, the craniodorsal section of the mesh will be fixed to the cervix. After prior hysterectomy or simultaneous hysterectomy performed with this surgery, the mesh will be attached to both sides of the vaginal apex.

Step 07: suture

- Last, the vaginal wound will be fully closed and sutured using absorbable sutures.
- Tamponade the vagina after the surgery.

Contraindications


- TiLOOP® meshes alter the expansion characteristics of surrounding tissue after implantation. Therefore, they should not be used in pregnant women or women whose family planning is not finished. Furthermore TiLOOP® meshes should not be applied on adolescent patients.
- TiLOOP® meshes should not be implanted in areas with active or latent infection.
- Do not use in the case of malignant situations of the vagina, the cervix, or the uterus.
- Patients at radiotherapy in the area of mesh implantation should not be treated with TiLOOP® meshes.
- Do not use on patients with intolerance against the product or known multiple or serious allergies.

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