
**Post-operative
wound drainage
Koelner Drainage**

The Koelner Drainage is a reservoir for post-operative wound drainage without suction. It allows the drainage of wound areas in which Redon drains are used without negative pressure.

www.pfmmedical.com



The **Koelner Drainage** is a reservoir for post-operative wound drainage without suction. It allows the drainage of wound areas in which Redon drains are used without negative pressure.

Detail View



Benefits

Hygienic 3 µm-outlet filter

When used as a gravity drainage the 3 µm-outlet filter allows easy escape of excess air in the reservoir.

OR-Version with universal Redon connector

The tube with transparent universal Redon connector CH 06 - CH 18 for connection to the drain is pre-attached.

Details

- ▶ Luer-Lock connection
- ▶ Crystal-clear material for reliable qualitative assessment of secretions
- ▶ Rough and fine graduation for at-a-glance information about current fill levels
- ▶ Quick-release fasteners included
- ▶ Double packed, sterile
- ▶ Redon bottle DEHP-free
- ▶ Latex-free

Technical Data

- ▶ **Bottle volume:** 400 ml

Ordering Information

REF	Version	Connection	PU
20384	Ward	Luer-Lock	40
20385	OP	Luer-Lock	40

Contact

Do you have any questions? Our Customer Solutions Team would be pleased to advise you.

service@pfmmedical.com
T +49 (0)2236 9641-220
F +49 (0)2236 9641-51

www.pfmmedical.com - Your source of information for pfm medical ag products

pfmmedical

pfm medical ag
Wankelstraße 60
50996 Köln, Germany

Certified according to
DIN EN ISO 13485